

So making CDS is a double infusion process and the recipe is this:

1. Put 300 milliliters of distilled water into the hermetic jar.
2. Place the shot glass into the jar.
3. Use the 10 ml syringe to measure 5 milliliters of Part A (22.4% sodium chlorite solution) and put it into the shot glass.
4. Use the 10 ml syringe to measure 5 milliliters of Part B (5% hydrochloric acid solution) and put it into the shot glass.
5. Close the lid down on the hermetic jar and latch it.
6. Put the jar away in a closet or little used cabinet for 24 hours.
7. 24 hours later, remove the jar from the cabinet, open the jar, remove the shot glass and close and latch the jar. Be careful to avoid the fumes.
8. Dump the contents of the shot glass into the garbage or store it in a plastic bottle for use as disinfectant.
9. Use the 10 ml syringe to measure 5 milliliters of Part A (22.4% sodium chlorite solution) and put it into the shot glass.
10. Use the 10 ml syringe to measure 5 milliliters of Part B (5% hydrochloric acid solution) and put it into the shot glass.
11. Re-open the hermetic jar and place the shot glass in it and close the lid and latch it.
12. Put the jar away for another 24 hours.